

Introduction to Project 1: Mapping relaxation times with MRI


Quantitative and Functional Imaging
BME 4420/7450
Fall 2022

Measuring relaxation times with MRI

- Map relaxation times using
 - Multiple echo time measurements of $T_2 = R_2^{-1}$
 - Inversion recovery measurements of $T_1 = R_1^{-1}$

Multiple echo time measurements of T_2


- Acquire a series of images, each at a different time after the tipping B_1 pulse
- Signal intensity decays exponentially with time constant T_2 .



Multiple echo time measurements of T_2


- Fit signal decay in each pixel to an exponential

$$S(T_E) = S_0 \cdot e^{-T_E/T_2}$$



Inversion recovery measurements of T_1


- Apply an 'inversion' B_1 pulse to orient M along $-Z$
- Wait some delay (inversion) time, T_I , during which M_z recovers toward M_0
- Tip spins into the transverse plane and measure magnetization immediately (indirect measure of M_z)



Inversion recovery measurements of T_1

- Acquire a series of images each with a different inversion time
- For each pixel location, fit signal recovery to estimate T_1

$$M_z(TI) = M_0 \cdot (1 - 2e^{-TI/T_1})$$



Goals of project

- Calculate and interpret T_2 and T_1 maps
- Compare T_1 and T_2 for tissue classification
- Practice basic MATLAB skills for image analysis

Making a head mask

- Identify pixels in background of image
- Compare pixel intensities to a threshold at 10% of maximum:

```
mask_m = (image_m > 0.1*max(image_m(:)));  
figure  
imagesc(mask_m)  
colormap(gray)  
axis image  
axis off  
title('Head mask')
```

- mask_m has 1's in head

